

INSTRUCTIONS TO OBTAIN DNA SAMPLE

1. Open a package containing sterile cotton swab (use rubber gloves if possible, do not touch swab end).
2. Touching only the handle of the swab, have the child open their mouth. Carefully rub the swab end on the inside of the child's mouth, inside the cheek.
3. Place swab aside, on a something, so that the swab is not touching anything. Let the swab dry for 10 minutes.
4. Seal swab in a **paper** envelope. Use self sealing envelope or tape to close. **Do not lick envelope** to seal it.
5. Label and store the envelope in a secure place.